

Product Name: Norwegian Gold Super Critical Omega

Product Class: Ultra Refined

Product Type: EE Capsules

Product Availability: Canada

Product Summary: Norwegian Gold Super Critical Omega Batch 34854 is a clean and highly concentrated omega-3 fish oil product. Each 1160 mg sample contains a minimum 80% omega-3 concentration as combined EPA and DHA, as well as gelatin, glycerin, lipase and mixed tocopherols. This product receives a FIVE out of FIVE Star IFOS rating based on:



- 1 Star for Passing All CRN/WHO Testing Categories - YES
- 1 Star for Confirmation of Label Claim for All Active Ingredients - YES
- 1 Star for Oxidation Level Less Than 75% of CRN Standard - YES
- 1 Star for PCB Levels Less Than 50% of CRN Standard - YES
- 1 Star for Dioxin Levels Less Than 50% of WHO Standard - YES

Batch / Lot #: 34854

Expiry Date: 02/2013

Recommended Daily Allowance: 1 Capsules

Company Name: Renew Life Canada

Company Phone #: 1-800-485-0960

Company Website: [www.renewlife.ca](http://www.renewlife.ca)





# CONSUMER REPORT

### CATEGORY 1 CONCENTRATION - Omega -3 Product Content vs. Product Label Claim

EPA

DHA

TOTAL OMEGA-3

| Product Label Claim | Batch Results | IFOS Compliance |
|---------------------|---------------|-----------------|
| 645 mg/1000mg       | 659 mg/1000mg | YES             |
| 100 mg/1000mg       | 131 mg/1000mg | YES             |
| 840 mg/1000mg       | 866 mg/1000mg | YES             |

### CATEGORY 2 PURITY, SAFETY AND CLEANLINESS - PCB's, Mercury, Dioxins and Furans

PCB

Dioxins and Furans

Dioxin Like PCB's

| CRN/GOED Standard | Detection Limit | Batch Result | IFOS Compliance |
|-------------------|-----------------|--------------|-----------------|
| < 90 ppb          | < ppt           | 3 ppb        | YES             |
| WHO Standard      | Detection Limit | Batch Result | IFOS Compliance |
| < 2 ppt           | < ppt           | 0.5 ppt      | YES             |
| < 3 ppt           | < ppt           | 0.85 ppt     | YES             |

### CATEGORY 3 STABILITY - Peroxide, Anisidine, Total Oxidation and Acid Value

Peroxide

Anisidine

Total Oxidation

Acid Value

| CRN/GOED Standard | Batch Results  | IFOS Compliance |
|-------------------|----------------|-----------------|
| < 5 meq/kg        | 1.78 meq/kg    | YES             |
| < 20 meq/kg       | 14.70 meq/kg   | YES             |
| < 26 meq/kg       | 18.26 meq/kg   | YES             |
| < 3.0 mg KOH/g    | < 2.0 mg KOH/g | YES             |

### CATEGORY 4 HEAVY METALS and MERCURY Lead, In-organic Arsenic, Cadmium

Mercury (Hg)

Lead (Pb)

Arsenic (As)

Cadmium (Cd)

| CRN/GOED Standard | Batch Results | IFOS Compliance |
|-------------------|---------------|-----------------|
| 0.1 ppm           | < 0.01 ppm    | YES             |
| 0.1 ppm           | < 0.01 ppm    | YES             |
| 0.1 ppm           | < 0.1 ppm     | YES             |
| 0.1 ppm           | < 0.01 ppm    | YES             |

